

TGFBR3 Rabbit pAb

货号: **B14283**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	93kDa
实测分子量	110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y,Mouse brain,Mouse lung,Rat brain
细胞定位	Cell membrane,Secreted,Single-pass type I membrane protein,extracellular space
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 140-340 of human TGFBR3 (NP_003234.2).
序列	ANFSLTAETEERNFPHGNEHLLNWARKEYGAVTSFTELKIARNIYIKVGEDQVFPPKCNIGKNFLSLNYLAEYLQPKAAEGCVMSSQPQNEEVHIIELITPNSNPYSAFQVDITIDIRPSQEDLEVVKNLILILKCKKSVNWVIKSFQVKGSLKIIAPNSIGFGKESERSMTMTKSIRDDIPSTQGNLVKWALDNGYSPIT

靶点信息

研究背景	This locus encodes the transforming growth factor (TGF)-beta type III receptor. The encoded receptor is a membrane proteoglycan that often functions as a co-receptor with other TGF-beta receptor superfamily members. Ectodomain shedding produces soluble TGFBR3, which may inhibit TGFB signaling. Decreased expression of this receptor has been observed in various cancers. Alternatively spliced transcript variants encoding different isoforms have been identified for this gene.
------	--

基因ID	7049
基因名	TGFBR3
Swiss	Q03167
别名	TGFBR3;BGCAN;betaglycan

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn