

SLC9A3R1 Rabbit pAb

货号: B14307

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Rattus norvegicus
应用	WB
推荐浓度	WB: 1:1000 - 1:4000
理论分子量	22kDa/38kDa
实测分子量	50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,293T,HT-29,MCF7,Mouse liver,Mouse kidney,Rat liver
细胞定位	Apical cell membrane,Cell projection,Cytoplasm,Endomembrane system,Peripheral membrane protein,filopodium,microvillus,ruffle
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 189-358 of human SLC9A3R1 (NP_004243.1).
序列	EASGLRAQDRIVEVNGVCMEGKQHGDVVS AIRAGGDET KLLVVDRETDEFFKKCRVIP SQEHLNGPLPVPFTNGEIQKEN SREALAEAALES PRPALVRSASSDTSEELNSQDSPPKQDSTAPSSTSSSDPILDFNISLAMAKERAHQKRSSKRAPQMDW SKKNELFSNL

靶点信息

研究背景	This gene encodes a sodium/hydrogen exchanger regulatory cofactor. The protein interacts with and regulates various proteins including the cystic fibrosis transmembrane conductance regulator and G-protein coupled receptors such as the beta2-adrenergic receptor and the parathyroid hormone 1 receptor. The protein also interacts with proteins that function as linkers between integral membrane and cytoskeletal proteins. The protein localizes to actin-rich structures including membrane ruffles, microvilli, and filopodia. Mutations in this gene result in hypophosphatemic nephrolithiasis/osteoporosis type 2, and loss of heterozygosity of this gene is implicated in breast cancer.
基因ID	9368
基因名	SLC9A3R1
Swiss	O14745
别名	SLC9A3R1;EBP50;NHERF;NHERF-1;NHERF1;NPHLOP2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn