

Phospho-FOXO3-S294 Rabbit pAb

货号: B14330

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Rattus norvegicus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	48kDa/71kDa
实测分子量	82kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	NIH/3T3
细胞定位	Cytoplasm,Nucleus,cytosol
纯化	Affinity purification

抗原信息


抗原信息	A synthetic phosphorylated peptide around S294 of human FOXO3 (NP_001446.1).
序列	PGSPT

靶点信息

研究背景	This gene belongs to the forkhead family of transcription factors which are characterized by a distinct fork head domain. This gene likely functions as a trigger for apoptosis through expression of genes necessary for cell death. Translocation of this gene with the MLL gene is associated with secondary acute leukemia. Alternatively spliced transcript variants encoding the same protein have been observed.
基因ID	2309
基因名	FOXO3

Swiss	O43524
别名	FOXO3;AF6q21;FKHRL1;FKHRL1P2;FOXO2;FOXO3A;FoxO3a

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn