

PPP1R15A Rabbit pAb

货号: B14332

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Danio rerio
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	49kDa/50kDa 56kDa/73kDa
实测分子量	110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,U-937,HeLa,HepG2,Jurkat,mouse spleen,mouse lung,mouse thymus
细胞定位	cytoplasm,cytosol,endoplasmic reticulum,Golgi apparatus,mitochondrial outer membrane,mitochondrion, protein phosphatase type 1 complex
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 375-674 of human PPP1 R15A (NP_055145.3).
序列	ASSSTPATGVFLKSWVYQPGEDTEEEEDSDTGSADEREAETSASTPPASAFLKAWVYRPGEDTEEEEDVDSEDK EDDSEAAALGEAESDPHPSPDQRAHFRGWGYRPGKETEEEEAAEDWGEAEPFVRVAIYVGEKPPPPWAPPRLPLRLQR RLKRPETPTHDPDPETPLKARKVRFSEKVTVHFLAVWAGPAQAARQGPWEQLARDRSRFARRITQAQEELSPCLTPAARA RAWARLRNPPLAPIALTQTLPSSSPVSPVQTTPLSQAVATPSRSSAAAAAALDLSGRRG

靶点信息

研究背景	This gene is a member of a group of genes whose transcript levels are increased following stressful growth arrest conditions and treatment with DNA-damaging agents. The induction of this gene by ionizing radiation occurs in certain cell lines regardless of p53 status, and its protein response is correlated with apoptosis following ionizing radiation.
基因ID	23645
基因名	PPP1R15A
Swiss	O75807
别名	PPP1R15A;GADD34

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn