

LARP6 Rabbit pAb

货号: B14340

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IF: Mus musculus WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	9kDa/54kDa
实测分子量	54kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Mouse brain,Rat brain
细胞定位	cytoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-320 of human LARP6 (NP_060827.2).
序列	MAQSGGEARPGPKTAVQIRVAIQEAEDVDELEDEEEGAETRGAGDPARYLSPGWGSASEEEPSRGHSGTTASGGENERE DLEQEWKPPDEELIKKLVDQIEFYFSDENLEKDAFLLLKHVRRNKLGYVSVKLLTSFKKVKHLTRDWRTTAHALKYSVVLELN EDHRKVRRTTPVPLFPNENLPSKMLLVYDLYLSPKLWALATPQKNGRVQEKVMEHLLKLFGTFGVISSVRILKPGRELPPDIR RISSRYSQVGTQECAIVEFEEVEAAIKAHEFMITESQGKENMKAVLIGMKPPKKKPAKDKNHDEEPTASIHNLKSL

靶点信息

研究背景	Regulates the coordinated translation of type I collagen alpha-1 and alpha-2 mRNAs, CO1A1 and CO1A2. Stabilizes mRNAs through high-affinity binding of a stem-loop structure in their 5' UTR. This regulation requires VIM and MYH10 filaments, and the helicase DHX9.
------	--

基因ID	55323
基因名	LARP6
Swiss	Q9BRS8
别名	LARP6;ACHN

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn