

ZMIZ2 Rabbit pAb

货号: B14350

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens IHC: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	115kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,Jurkat,HeLa,293T,mouse brain,rat brain
细胞定位	mitochondrion,nuclear replication fork,nucleoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human ZMIZ2 (NP_001287888).
序列	MNSMNPMPALPPAPHGDGSFAYESVPWQQSATQPAGSLSVTTVWGVGNATQSQCLGQQAF AEGGANKGYVQQG VYSRGGYPGAPGFTTGYAGGPGGLGLPSHAARPST

靶点信息

研究背景	ZMIZ2 and ZMIZ1 (MIM 607159) are members of a PIAS (see MIM 603566)-like family of proteins that interact with nuclear hormone receptors. ZMIZ2 interacts with androgen receptor (AR; MIM 313700) and enhances AR-mediated transcription (Huang et al., 2005 [PubMed 16051670]).
基因ID	83637

基因名	ZMIZ2
Swiss	Q8NF64
别名	ZMIZ2;NET27;TRAFIP20;ZIMP7;hZIMP7

产品验证



Western blot analysis of ZMIZ2 expressed in mouse brain, rat brain, U-87MG using ZMIZ2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn