

TET2 Rabbit pAb

货号: B14351

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:100 - 1:500
理论分子量	130kDa/133kDa/223kDa
实测分子量	280KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y,293T
细胞定位	nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of TET2 (NP_001120680.1).
序列	MEQDRTNHVEGNRLSPFLIPSPPICQTEPLATKLQNGSPLPERAHPEVNGDTKWHSFKSYYGIPCMKGSQNSRVSPDFTQESRGYSKCLQNGGIKRTVSEPSLSGLLQIKLKQDQKANGERRNFGV SQERNPGESSQPNVSDLSDDKESVSSVAQENAVKDFTSFSTHNCSGPNPELQILNEQEGKSANYHDKNIVLLKNKAVLMPNGATVSASSVEHTHGELLEKTL SQYYPDCVSI AVQKTTSHINAINSQATNELSCEITHPSHTSGQINSAQTSNSELPPKPAAVVSEACDADDADNASKLAAMLNTCSFQKPEQLQQKSVFEICPSAENNIQGTTKLASGE

靶点信息

研究背景	The protein encoded by this gene is a methylcytosine dioxygenase that catalyzes the conversion of methylcytosine to 5-hydroxymethylcytosine. The encoded protein is involved in myelopoiesis, and defects in this gene have been associated with several myeloproliferative disorders. Two variants encoding different isoforms have been found for this gene.
------	--

基因ID	54790
基因名	TET2
Swiss	Q6N021
别名	MDS;KIAA1546;TET2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn