

# RPF2 Rabbit pAb

货号: B14389

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                      |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  | <b>WB:</b> Homo sapiens                                                                              |
| 应用    | <a href="#">WB</a>                                                                                   |
| 推荐浓度  | <b>WB:</b> 1:100 - 1:500                                                                             |
| 理论分子量 | 35kDa                                                                                                |
| 实测分子量 | 36KDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | HeLa                                                                                                 |
| 细胞定位  | nucleolus,nucleoplasm                                                                                |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                                                                                                                                                                                                                                                |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-306 of human RPF2 (N P_115570.1).                                                                                                                                                                                                              |
| 序列   | MDTLDRVVKPKTKRAKRFLEKREPKLNENIKNAMLIKGGNANATVTKVLKDVYALKKPYGVLYKKKNITRPFEDQTSLEFFS<br>KKSDCSLFMFGSHNKKRPNNLVIGRMYDYHVLDMIELGIENFVSLKDIKNSKCPEGTKPMLIFAGDDFDVTEDYRRLKSLLI<br>DFFRGPTVSNIRLAGLEYVLHFTALNGKIYFRSYKLLKKSGCRTPRIELEEMGPSLDLVLRRTHLASDDLYKL SMKMPKALK<br>PKKKKNISHDTFGTTYGRIHMQKQDLSKLQTRKMKGLKKRPAERITEDHEKKSKRIKKN |

靶点信息

|      |                                                                                                                                           |
|------|-------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | Involved in ribosomal large subunit assembly. May regulate the localization of the 5S RNP/5S ribonucleoprotein particle to the nucleolus. |
| 基因ID | 84154                                                                                                                                     |

|       |                      |
|-------|----------------------|
| 基因名   | RPF2                 |
| Swiss | Q9H7B2               |
| 别名    | RPF2;BXDC1;bA397G5.4 |

产品验证



Western blot analysis of RPF2 expressed in HeLa using RPF2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)