

CIZ1 Rabbit pAb

货号: B14396

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens IF: Homo sapiens IHC: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	100kDa
实测分子量	100kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	mouse testis
细胞定位	nucleoplasm,nucleus,plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 250-440 of human CIZ1 (NP_001124490).
序列	PGQLQVKAQPQARMTVPKQTQTPDLLPEALEAQVLPRFQPRVLQVQAQVQSQTQPRIPSTDTQVQPKLQKQAQTQTSPEHLVLQQKQVQPQLQQAEPQKQVQPQVHTQAQPSVQPQEHPAQVSVQPPEQTHEQPHTQPQVSLLAPEQTPVVVHVCGLEMPDDAVEAGGGMEKTLPEPVGTQVSMEEIQN

靶点信息

研究背景	The protein encoded by this gene is a zinc finger DNA binding protein that interacts with CIP1, part of a c omplex with cyclin E. The encoded protein may regulate the cellular localization of CIP1. Several transcrip t variants encoding different isoforms have been found for this gene.
------	---

基因ID	25792
基因名	CIZ1
Swiss	Q9ULV3
别名	CIZ1;LSFR1;NP94;ZNF356

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn