

NTNG1 Rabbit pAb

货号: B14408

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: mice
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	60kDa
实测分子量	60kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,mouse brain,mouse lung,mouse uterus,rat brain
细胞定位	extracellular region,glutamatergic synapse,laminin complex,plasma membrane,presynaptic active zone membrane,Schaffer collateral - CA1 synapse
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 29-300 of human NTNG 1 (NP_055732.2).
序列	HYDLCKTQIYTEEGKVWDYMACQPESTDMTKYLKVKLDPPDITCGDPPETFCAMGNPYMCNNECDASTPELAHPPELMFDFEGRHPSTFWQSATWKEYPKPLQVNITLSWSKTIELTDNIVITFESGRPDQMILEKSLDYGRTWQPYQYYATDCLDAFHMDPKSVKDLSQHTVLEIICTEEYSTGYTTNSKIIHFEIKDRFAFFAGPRLRNMASLYGQLDTTKKLRDFFTVDLRI RLLRPVAGEIFVDELHLARYFYAISDIKVRGRCKCN

靶点信息

研究背景	This gene encodes a preproprotein that is processed into a secreted protein containing eukaroytic growth factor (EGF)-like domains. This protein acts to guide axon growth during neuronal development. Polymorphisms in this gene may be associated with schizophrenia. Alternative splicing results in multiple transcript variants encoding distinct isoforms.
------	---

基因ID	22854
基因名	NTNG1
Swiss	Q9Y2I2
别名	Lmnt1;NTNG1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn