

POU3F4 Rabbit pAb

货号: B14414

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Crassostrea gigas
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	42kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,Mouse brain,Rat brain
细胞定位	nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human POU3F4 (NP_000298.3).
序列	MATAASNPYSILSSTSLVHADSAGMQQGSPFRNPQKLLQSDYLQGVPSNGHPLGHHWVTSLSDGGPWSSTLATSPLDQ QDVKPGREDLQLGAIHHRSPhVAHHSPHTNHPNAWGASPAPNPSITSSGQPLNVYSQPGFTVSGMLEHGGLTPPPAAA SAQSLHPVLRPPDHGELGSHHCQDHSDEETPT

靶点信息

研究背景	This gene encodes a member of the POU-III class of neural transcription factors. This family member plays a role in inner ear development. The protein is thought to be involved in the mediation of epigenetic signals which induce striatal neuron-precursor differentiation. Mutations in this gene are associated with X chromosome-linked nonsyndromic mixed deafness. [provided by RefSeq, Dec 2012]
基因ID	5456

基因名	POU3F4
Swiss	P49335
别名	BRAIN-4;BRN-4;BRN4;DFN3;DFNX2;OCT-9;OTF-9;OTF9;POU3F4

产品验证



Western blot analysis of POU3F4 expressed in U-87MG,Mouse brain,Rat brain using POU3F4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn