

SP140 Rabbit pAb

货号: B14426

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	98kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse spleen,Mouse thymus,Rat thymus
细胞定位	fibrillar center,mitochondrion,nucleus,PML body
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human SP140 (NP_001265380).
序列	MAQQGQQGQMASGDSNLNFRMVAEIQNVEGQNLQEVCPEIFRFFRENKVEIASAITRPFPFLMGLRDRSFISEQMYEH FQEAFRNLVPVTRVMYCVLSELEKTFGWSHLEALFSRINLMAYPDLNEIYRSFQNVCEHSPLQMNNVNDLEDPRLLPYG KQENSNACHEMDDIAPVQE

靶点信息

研究背景	This gene encodes a member of the SP100 family of proteins, which are share common domains including an N-terminal homogeneously staining region domain followed by a SP100/autoimmune regulator/NucP41/P75/deformed epidermal autoregulatory factor domain, a plant homeobox zinc finger, and a bromodomain. The encoded protein is interferon-inducible and is expressed at high levels in the nuclei of leukocytes. Variants of this gene have been associated with multiple sclerosis, Crohn's disease, and chronic lymphocytic leukemia. Alternative splicing results in multiple variants. [provided by RefSeq, Aug 2016]
------	---

基因ID	11262
基因名	SP140
Swiss	Q13342
别名	LYSP100;LYSP100-A;LYSP100-B;SP140

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn