

THSD4 Rabbit pAb

货号: B14435

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	112kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y
细胞定位	extracellular exosome
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 500-700 of human THSD4 (NP_079093.2).
序列	AEGPTNEILDVYMIHQPNPGVHYEYVIMGTNAISPQVPPHRRPGEPFNGQMVTTEGRSQEEGEQKGRNEEKEDLRGEAPE MFTSESAQTFPVRHPDRFSPHRPDNLVPPAPQPPRRSRDHNWKQLGTTECSTTCGKGSQYPIFRCVHRSTHEEAPESYCD SSMKPTPEEEPCNIFPCPAFWDIGEWSECSKTCGLGMQHRQ

靶点信息

研究背景	Promotes FBN1 matrix assembly. Attenuates TGFB signaling, possibly by accelerating the sequestration of large latent complexes of TGFB or active TGFB by FBN1 microfibril assembly, thereby negatively regulating the expression of TGFB regulatory targets, such as POSTN (By similarity).
基因ID	79875

基因名	THSD4
Swiss	Q6ZMP0
别名	ADAMTSL-6;ADAMTSL6;FVSY9334;PRO34005;THSD4

产品验证



Western blot analysis of THSD4 expressed in SH-SY5Y using THSD4 Rabbit pAb at 1:1000. Secondary anti body: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% no n-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn