

GPR35 Rabbit pAb

货号: B14440

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	34kDa
实测分子量	34kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,mouse stomach,mouse liver
细胞定位	plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human GPR35 (NP_00529 2.2).
序列	QQWTETRIYMTNLAVADLCLLCTLPFVLHSLRDTSDTPLCQLSQGIYLTNRYMSISLVTAIAVDRYVAVRHPLRARGLRSPR QAAAVCAVLWVLVIGSLVA

靶点信息

研究背景	GPR35 (G protein-coupled receptor 35) is a protein-coding gene. Diseases associated with GPR35 include primary sclerosing cholangitis. GO annotations related to this gene include G-protein coupled receptor activity. An important paralog of this gene is GPR18.
基因ID	2859

基因名	GPR35
Swiss	Q9HC97
别名	GPR35

产品验证



Western blot analysis of GPR35 expressed in Mouse brain,mouse stomach,mouse liver using GPR35 Rab bit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30 ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn