

XRN2 Rabbit pAb

货号: B14510

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	IF: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	108kDa
实测分子量	108kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	PC-3,A-549
细胞定位	nucleolus,nucleoplasm,nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-950 of human XRN2 (Q9H0D6).
序列	MGVPAFFRWLSRKYPsiivncveekpkECNGVKIPVDASKPNPNdVEFDNLYLDMNGIIHPCTHPEDKPAPKNEDEMMVAIF EYIDRLFSIVRPRRLLYMAIDGVAPRAKMNQQRsRRFRASKEGMEAAVEKQRVREEILAKGGFLPPEEIKERFDSNCITPGTE FMDNLAKCLRYIIADRLNNDPGWKNLTVILSDASAPGEGEHKIMDYIRRQRAQPNHDPNTHHCLCGADADLIMLGLATHE PNFTIIEEEFKPNKPKPCGLCNQFGHEVKDCEGLPREKKGKHDELADSLPCAEGEFIFLRLNVLREYLERELTMASLPFTFDV ERSIDDWVFMCFVGNDFLPHLPsLEIRENAIDRLVNIYKNVVHKTGGYLTESGYVNLQRVQMIMLAVGEVEDSIFKKRKD DEDSFRRRQKEKRKRMRKDQPAFTPSGILTPHALGSRNSPGSQVASNPRQAAYEMRMQNNSSPSISPNTSFTSDGSPSP LGGIKRKAEDSDSEPEPEDNVRLWEAGWKQRYYKNKFDVDAADEKFRRKVVQSYVEGLCWVLRYYYYQGcASWKWYYp FHYAPFASDFEGiADMPsDFEKGTKPKFPLEQLMGVFPAASGNFLPPSWRKLMsDPDSSiIDFYPEDFAIDLNGKKYAWQ GVALLPFVDERRLRaALEEVYPDLTPETRRNSLGGDVLfVGKHHPLHDFILELYQTGSTEPVEVPPELCHGIQgKfSLDEE AILPDQIVCSPVPMRLDLTQNTVVSINFKDpQFAEDYIFKAVMLPGARKPAAVLKPSDWEKSSNGRQWKpQLGFNRDRRP VHLDQAAfRTLGHVMPRGSGTGIYSNAAPPpVTYQGNLYRPLLRGQAQIPKLMSNMRPQDSWRGPPpLFQQQRfDRGV GAEPllPWNRLQTQNAAFQPNQYQMLAGPGGYPPRRDDRGRQGYPREGRKYPLPPPSGRYNWN

靶点信息

研究背景	This gene encodes a 5'-3' exonuclease that promotes transcription termination at cotranscriptional cleavage sites. Alternative splicing results in multiple transcript variants encoding different isoforms.
基因ID	22803
基因名	XRN2
Swiss	Q9H0D6
别名	XRN2

产品验证



Western blot analysis of XRN2 expressed in PC-3,A-549 using XRN2 Rabbit pAb at 1:1000. Secondary anti body: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% no n-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn