

SH2D2A Rabbit pAb

货号: B14521

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	RT-PCR: Mus musculus WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	39kDa/40kDa/42kDa/44kDa
实测分子量	41kDa/43kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,K-562,Mouse spleen,Rat thymus
细胞定位	cytoplasm,cytosol
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 250-350 of human SH2D2A (NP_003966.2).
序列	PAKPQLPPEVYTIPVPRHRPAPRPKPSNPIYNEPDEPIAFYAMGRGSPGEAPSNIYVEVEDEGLPATLGHPVLRKSWSRPVPG GQNTGGSQ LHSENSVIGQ

靶点信息

研究背景	This gene encodes an adaptor protein thought to function in T-cell signal transduction. A related protein in mouse is responsible for the activation of lymphocyte-specific protein-tyrosine kinase and functions in downstream signaling. Alternative splicing results in multiple transcript variants.
基因ID	9047

基因名	SH2D2A
Swiss	Q9NP31
别名	F2771;SCAP;TSAD;VRAP;SH2D2A

产品验证



Western blot analysis of SH2D2A expressed in Mouse spleen,Rat thymus,U-87MG using SH2D2A Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 12 0s.

实验步骤

访问官网浏览详情: www.ablybio.cn