

PNPO Rabbit pAb

货号: B14522

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	7kDa/25kDa/27kDa/29kDa
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,293T,Mouse brain,Mouse kidney,Mouse liver,Rat brain,Rat kidney,Rat liver
细胞定位	cytosol
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-261 of human PNPO (NP_060599.1).
序列	MTCWLRGVTATFGRPAEWPGYLSHLCGRSAAMDLGPMRKSYPGRDREAFEETHLTSLDPVKQFAAWFEEAVQCPDIGEANAMCLATCTRDGKPSARMLLLKGFGKDGFRFFTNFESRKGKELDSNPFASLVFYWEPLNRQVRVEGPVKKLPAAAAECYFHSRPAKSSQIGAVVSHQSSVIPDREYLRKKNEELEQLYQDQEVKPKSWGGYVLYPQVMEFWQGGQTNRLLHDRVFRRLPTGDSPLGPMTHRGEEDWLYERLAP

靶点信息

研究背景	The enzyme encoded by this gene catalyzes the terminal, rate-limiting step in the synthesis of pyridoxal 5'-phosphate, also known as vitamin B6. Vitamin B6 is a required co-factor for enzymes involved in both homocysteine metabolism and synthesis of neurotransmitters such as catecholamine. Mutations in this gene result in pyridoxamine 5'-phosphate oxidase (PNPO) deficiency, a form of neonatal epileptic encephalopathy.
基因ID	55163
基因名	PNPO
Swiss	Q9NVS9
别名	HEL-S-302;PDXPO;PNPO

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn