

VAMP7 Rabbit pAb

货号: B14693

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IF/ICC
推荐浓度	WB: 1:100 - 1:500 IF/ICC: 1:50 - 1:200
理论分子量	20kDa/24kDa/30kDa
实测分子量	25KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-1080,Mouse brain,Rat brain
细胞定位	cytoplasm,extracellular exosome,neuron projection,plasma membrane,SNARE complex,synapse,trans-Golgi network
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 11-70 of human VAMP7 (P51809).
序列	GTTILAKHAWCGGNFLEVTEQILAKIPSENNKLTYSHGNYLFHYICQDRIVYLCITDDDF

靶点信息

研究背景	This gene encodes a transmembrane protein that is a member of the soluble N-ethylmaleimide-sensitive f actor attachment protein receptor (SNARE) family. The encoded protein localizes to late endosomes and lysosomes and is involved in the fusion of transport vesicles to their target membranes. Alternate splicing results in multiple transcript variants.
基因ID	6845

基因名	VAMP7
Swiss	P51809
别名	SYBL1; TIVAMP; VAMP-7; TI-VAMP

产品验证



Western blot analysis of VAMP7 expressed in HT-1080,Mouse brain,Rat brain using VAMP7 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn