

RAB24 Rabbit pAb

货号: B14918

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Human liver cancer
应用	WB IP
推荐浓度	WB: 1:500 - 1:2000 IP: 1:50 - 1:100
理论分子量	23kDa
实测分子量	23kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-1080,MCF7,SH-SY5Y,293T,Mouse brain
细胞定位	Cytoplasm,Lipid-anchor,Membrane,cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-203 of human RAB24 (NP_001026847.1).
序列	MSGQRVDVKVVMVGKEYVGKTSLVERYVHDRFLVGPYQNTIGAAFVAKVMSVGDRTVTLGIWDTAGSERYEAMSRIYYR GAKAAIVCYDLTDSSSFERAKFWVKELRSLEEGCQIYLCGTKSDLLEEDRRRRRVDFHDVQDYADNIKAQLFETSSKTGQ SVDELFGKVAEDYVSAAFQVMTEDKGVDLGQKPNPYFYSCCHH

靶点信息

研究背景	RAB24 is a small GTPase of the Rab subfamily of Ras-related proteins that regulate intracellular protein trafficking (Olikkonen et al., 1993 [PubMed 8126105]).[supplied by OMIM, Aug 2009]
基因ID	53917

基因名	RAB24
Swiss	Q969Q5
别名	RAB24

产品验证



Western blot analysis of RAB24 expressed in Mouse brain, HT-1080, MCF7 using RAB24 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn