

# DTX2 Rabbit pAb

货号: B14922

产品信息

|       |                                                                                                        |
|-------|--------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                        |
| 宿主    | Rabbit                                                                                                 |
| 克隆性   | Polyclonal                                                                                             |
| 预测反应  | <b>WB:</b> 293T                                                                                        |
| 应用    | <a href="#">WB</a> <a href="#">IHC</a>                                                                 |
| 推荐浓度  | <b>WB:</b> 1:100 - 1:500<br><b>IHC:</b> 1:50 - 1:200                                                   |
| 理论分子量 | 62kDa/67kDa                                                                                            |
| 实测分子量 | 67KDa                                                                                                  |
| 形式    | Liquid                                                                                                 |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                           |
| 阳性对照  | A-549                                                                                                  |
| 细胞定位  | Cytoplasm,Nucleus                                                                                      |
| 纯化    | Affinity purification                                                                                  |

抗原信息


|      |                                                                                                                                                                                                                                                                                                     |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human DTX2 ( NP_065943.2).                                                                                                                                                                                   |
| 序列   | MAMAPSPSLVQVYTSPAAVAVWEWQDGLGTWHPYSATVCSFIEQQFVQQKGQRFLGSLAHSIPLGQADPSLAPYIIDL<br>PSWTQFRQDTGTMRAVRRHLFPQHSAPGRGVVWEWLSDDGSWTAYEASVCDYLEQQVARGNQLVDLAPLGYNNTVN<br>YTTHTQTNKTSSFCRSVRRQAGPPYPVTIIAPPGHTGVACSCHQCLSGSRTGPVSGRGRHSMTNLPAYPVPQHPPHRTAS<br>VFGTHQAFAPYINKPSLSGARSAPRLNTTNAWGAAPPSLGSQPLY |

靶点信息

|      |                                                                                                                 |
|------|-----------------------------------------------------------------------------------------------------------------|
| 研究背景 | DTX2 functions as an E3 ubiquitin ligase (Takeyama et al., 2003 [PubMed 12670957]).[supplied by OMIM, Nov 2009] |
| 基因ID | 113878                                                                                                          |

|       |            |
|-------|------------|
| 基因名   | DTX2       |
| Swiss | Q86UW9     |
| 别名    | DTX2;RNF58 |

产品验证



Western blot analysis of DTX2 expressed in A-549 using DTX2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)