

SNRPD2 Rabbit pAb

货号: B14923

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: 293T
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	12kDa/13kDa
实测分子量	13kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,HL-60,HepG2,HeLa,Jurkat,Rat testis,Rat heart
细胞定位	Cytoplasm,Nucleus,cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-118 of human SNRPD 2 (NP_004588.1).
序列	MSLLNPKPKSEMTPEELQKREEEFNTGPLSVLTQSVKNNTQVLINCRNNKLLGRVKAFDRHCNMVLENVKEMWTEVPKS GKGKKKSKPVNKDRYISKMFLRGDSVIVVLRNPLIAGK

靶点信息

研究背景	The protein encoded by this gene belongs to the small nuclear ribonucleoprotein core protein family. It is required for pre-mRNA splicing and small nuclear ribonucleoprotein biogenesis. Multiple transcript variant s encoding different isoforms have been found for this gene.
基因ID	6633

基因名	SNRPD2
Swiss	P62316
别名	SNRPD2;SMD2;SNRPD1;Sm-D2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn