

SPTLC1 Rabbit pAb

货号: B14924

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: HeLa
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	16kDa/52kDa
实测分子量	53kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SKOV3,MCF7,HepG2,Mouse kidney,Mouse lung,Mouse liver,Mouse thymus,Rat lung,Rat kidney
细胞定位	Endoplasmic reticulum membrane,Single-pass membrane protein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 37-310 of human SPTLC1 (NP_006406.1).
序列	RLLFSKTYKLQERSDLTVKEKEELIEEWQPEPLVPPVPKDHPALNYNIVSGPPSHKTVVNGKECINFASFNFLGLLDNPRVKA AALASLKKYGVGTCGPRGFYGTFDVHLDLEDRLAKFMKTEEAIISYGFATIASAIPAYSKRGDIVFDRAACFAIQKGLQAS RSDIKLFKHNDMADLERLLKEQIEDQKNPRKARVTRRFIVVEGLYMNTGTICPLPELVKLKYKYKARIFLEESLSFGVLGEH GRGVTEHYGINIDDIDLISANMENA

靶点信息

研究背景	This gene encodes a member of the class-II pyridoxal-phosphate-dependent aminotransferase family. The encoded protein is the long chain base subunit 1 of serine palmitoyltransferase. Serine palmitoyltransferase converts L-serine and palmitoyl-CoA to 3-oxosphinganine with pyridoxal 5'-phosphate and is the key enzyme in sphingolipid biosynthesis. Mutations in this gene were identified in patients with hereditary sensory neuropathy type 1. Alternatively spliced variants encoding different isoforms have been identified. Pseudogenes of this gene have been defined on chromosomes 1, 6, 10, and 13.
基因ID	10558
基因名	SPTLC1
Swiss	O15269
别名	SPTLC1;HSAN1;HSN1;LBC1;LCB1;SPT1;SPTI

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn