

# NLRX1 Rabbit pAb

货号: B14944

## 产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                      |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB IP                                                                                                |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IP:</b> 1:500 - 1:1000                                               |
| 理论分子量 | 101kDa/107kDa                                                                                        |
| 实测分子量 | 108kDa                                                                                               |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | MCF7,Mouse heart,Mouse kidney,Mouse liver,Mouse lung,Rat liver                                       |
| 细胞定位  | Mitochondrion outer membrane                                                                         |
| 纯化    | Affinity purification                                                                                |

## 抗原信息


|      |                                                                                                                                                                                                                                                                                               |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 87-360 of human NLRX1 (NP_001269073.1).                                                                                                                                                                         |
| 序列   | AEWFSRLPREERQFGPTFALDTVHVDPVIRESTPDELLRPPAELALEHQPPQAGLPPLALSQLFNPDACGRRVQTVVLYGTV<br>GTGKSTLVRKMVLDWCYGRLPAFELLIPFSCEDLSSLGPAPASLCQLVAQRYTPLKEVLPMAAAGSHLLFVLHGLEHLNLD<br>FRLAGTGLCSDPEEPQEPAAIIVNLLRKYMPLQASILVTTRPSAIGRIPSKYVGRYGEICGFSDTNLQKLYFQLRLNQPYCGYA<br>VGGSGVSATPAQRDHLVQMLSRNLEG |

## 靶点信息

|      |                                                                                                                                                                                                                                                                               |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | The protein encoded by this gene is a member of the NLR family and localizes to the outer mitochondrial membrane. The encoded protein is a regulator of mitochondrial antiviral responses. Three transcript variants encoding the same protein have been found for this gene. |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                                      |
|-------|--------------------------------------|
| 基因ID  | 79671                                |
| 基因名   | NLRX1                                |
| Swiss | Q86UT6                               |
| 别名    | NLRX1;CLR11.3;DLNB26;NOD26;NOD5;NOD9 |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)