

CSRP1 Rabbit pAb

货号: B14977

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: A172
应用	WB IP RIP
推荐浓度	WB: 1:500 - 1:2000 IP: 1:50 - 1:200 RIP: 1:20 - 1:50
理论分子量	20kDa
实测分子量	21kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,BT-474,HeLa,Mouse lung,Mouse brain
细胞定位	Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-193 of human CSRP1 (NP_004069.1).
序列	MPNWGGGKKCGVCQKTVYFAEEVQCEGNSFHKSCFLCMVCKKNLDSTTVAVHGEEIYCKSCYGKKYGPKGYGYGQGA GTLSTDKGESLGIKHEEAPGHRPTTNPNASKFAQKIGGSERCPRCSQAVYAAEKVIGAGKSWHKACFRCAKCGKGLESTT LADKDGEIYCKGCYAKNFGPKGFGFGQGAGALVHSE

靶点信息

研究背景	This gene encodes a member of the cysteine-rich protein (CSRP) family. This gene family includes a group of LIM domain proteins, which may be involved in regulatory processes important for development and cellular differentiation. The LIM/double zinc-finger motif found in this gene product occurs in proteins with critical functions in gene regulation, cell growth, and somatic differentiation. Alternatively spliced transcript variants have been described.
基因ID	1465
基因名	CSRP1
Swiss	P21291
别名	CSRP1;CRP;CRP1;CSRP;CYRP;D1S181E;HEL-141;HEL-S-286

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn