

DUSP15 Rabbit pAb

货号: B15056

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:200
理论分子量	14kDa/25kDa/26kDa/31kDa
实测分子量	32kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	PC-3,A-375,Mouse kidney,Rat kidney
细胞定位	cytosol,plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 10-90 of human DUSP15 (NP_542178.2).
序列	PGLYLGNFIDAKDLDDLGRNKITHIISIHESQPQLLDITYLRIPVADTPEVPIKKHFKECINFIHCCRLNGGNCLVHCFA

靶点信息

研究背景	The protein encoded by this gene has both protein-tyrosine phosphatase activity and serine/threonine-specific phosphatase activity, and therefore is known as a dual specificity phosphatase. This protein may function in the differentiation of oligodendrocytes. Alternative splicing results in multiple transcript variants.
基因ID	128853

基因名	DUSP15
Swiss	Q9H1R2
别名	DUSP15;C20orf57;VHY

产品验证



Western blot analysis of DUSP15 expressed in Mouse kidney,Rat kidney,PC-3 using DUSP15 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn