

FAM160B2 Rabbit pAb

货号: B15060

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC IP
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:500 - 1:1000
理论分子量	
实测分子量	82KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat testis
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 74-303 of mouse FAM16 OB2 (NP_919326.1).
序列	GEAGPCLEYLLQHKILETLCTLGKAEYPPGMRQQVFQFFSKVLSQVQHPLLHYLSVHRPVQKLLRLGGTVPGSLTEKEEVQ FTSVLCSKIQQDPELLAYILEGKKIIGKKKTARESTAPPKDIAGYRDKDCPHSDALNRDPGLDKEHCGVPALSIHLPAETEGP ENGPGESNLITSLGLCKSKKSRLALKAQENILLVSVASPAATYLTQSTSCCMAIAEHLQCQLYR

靶点信息

研究背景	Able to activate MAPK/ERK and TGFB signaling pathways (By similarity. May regulate the activity of genes involved in intestinal barrier function and immunoprotective inflammation. May play a role in cell proliferation (By similarity.
------	---

基因ID	239170
基因名	
Swiss	Q80YR2
别名	G430067P06Rik; Rai16

产品验证



Western blot analysis of FAM160B2 expressed in Mouse brain,Rat testis using FAM160B2 Rabbit pAb at 1 :1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn