

Phospho-HDAC8-S39 Rabbit pAb

货号: B15083

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IHC IF/ICC
推荐浓度	IHC: 1:50 - 1:100 IF/ICC: 1:100 - 1:200
理论分子量	15kDa/16kDa/17kDa/28kDa/29kDa/31kDa/41kDa
实测分子量	15kDa/16kDa/17kDa/28kDa/29kDa/31kDa/41kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息

抗原信息	A phospho specific peptide corresponding to residues surrounding S39 of human HDAC8
序列	

靶点信息

研究背景	Histones play a critical role in transcriptional regulation, cell cycle progression, and developmental events . Histone acetylation/deacetylation alters chromosome structure and affects transcription factor access to DNA. The protein encoded by this gene belongs to class I of the histone deacetylase family. It catalyzes t he deacetylation of lysine residues in the histone N-terminal tails and represses transcription in large mul ti-protein complexes with transcriptional co-repressors. Multiple transcript variants encoding different isof orms have been found for this gene.
基因ID	55869

基因名	HDAC8
Swiss	Q9BY41
别名	HDAC8;CDA07;CDLS5;HD8;HDACL1;MRXS6;RPD3;WTS

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn