

VPS51 Rabbit pAb

货号: B15109

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	72kDa/86kDa
实测分子量	100kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549
细胞定位	Golgi apparatus,Recycling endosome,trans-Golgi network
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human VPS51 (NP_037397.2).
序列	MAAAAAAGSPGSGPGDSPEGPEGEAPERRRKAHGMLKLYGLSEGEAAGRPAGPDPLDPTDLNGAHFDPEVYLDKLRRECP LAQLMDSETDMVRQIRALDSMQTLVYENYNKFISATDTIRKMKNDFRKMEDEMDRLATNMAVITDFSARISATLQDRHERIT KL AGVHALLRKLQFLFELPSRLTKCVELGAYGQAVRYQGRAQAVLQYQHLPSFRAIQDDCQVITARLAQQLRQRFREGGSGA PEQAECVELLLA

靶点信息

研究背景	This gene encodes a member of the vacuolar protein sorting-associated protein 51 family. The encoded p rotein is a component of the Golgi-associated retrograde protein complex which acts as a tethering factor for carriers in retrograde transport from the early and late endosomes to the trans-Golgi network. Alterna tive splicing results in multiple transcript variants.
------	---

基因ID	738
基因名	VPS51
Swiss	Q9UID3
别名	VPS51;ANG2;ANG3;C11orf2;C11orf3;FFR

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn