

SLC18A3 Rabbit pAb

货号: B15121

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IHC IF/ICC
推荐浓度	IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	56kDa
实测分子量	65KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 473-533 of human SLC18A3 (NP_003046.2).
序列	GLLTRSRSERDVLLDEPPQGLYDAVRLRERPVSQGDEPRSPPGPFDACEDDYNYYYTRS

靶点信息

研究背景	This gene is a member of the vesicular amine transporter family. The encoded transmembrane protein transports acetylcholine into secretory vesicles for release into the extracellular space. Acetylcholine transport utilizes a proton gradient established by a vacuolar ATPase. This gene is located within the first intron of the choline acetyltransferase gene.
基因ID	6572

基因名	SLC18A3
Swiss	Q16572
别名	SLC18A3;CMS21;VACHT

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn