

GTPBP8 Rabbit pAb

货号: B15125

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	13kDa/28kDa/32kDa
实测分子量	32kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Rat liver
细胞定位	mitochondrion
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human GTPBP8 (NP_054889.2).
序列	MAAPGLRLGAGRLFEMPAVLERLSRYNSTSQAFAEVLRLPKQQLRKLLYPLQEVERFLAPYGRQDLHLRIFDPSPEDIARAD NIFTATERNRIDYVSSAVRIDHAPDLRPEVCFIGRSNVGKSSLIKALFSLAPEVEVRVSKKPGHTKKMNFFKVGKHFTVVD MPGYGF

靶点信息

研究背景	Predicted to enable GTP binding activity and metal ion binding activity. Located in mitochondrion.
基因ID	29083

基因名	GTPBP8
Swiss	Q8N3Z3
别名	GTPBP8;HSPC135

产品验证



Western blot analysis of GTPBP8 expressed in Mouse liver,Rat liver using GTPBP8 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn