

D2HGDH Rabbit pAb

货号: B15128

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	32kDa/34kDa/56kDa
实测分子量	56kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Rat heart
细胞定位	mitochondrial matrix,mitochondrion
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 322-521 of human D2H GDH (NP_689996.4).
序列	GRHLHLASPVQESPFYVLIETSGSNAGHDAEKLGHFLEHALGSGLVTDGTMATDQRKVKMLWALRERITEALSRDGYVYK YDLSLPVERLYDIVTDLRLARLGPHAKHVVGYGHLGDGNLHLNVTAEAFSPSLAALPHVYEWTAGQQGSVSAEHGVGFR KRDVLGYSKPPGALQLMQQLKALLDPKGILNPYKTLPSQA

靶点信息

研究背景	This gene encodes D-2hydroxyglutarate dehydrogenase, a mitochondrial enzyme belonging to the FAD-binding oxidoreductase/transferase type 4 family. This enzyme, which is most active in liver and kidney but also active in heart and brain, converts D-2-hydroxyglutarate to 2-ketoglutarate. Mutations in this gene are present in D-2-hydroxyglutaric aciduria, a rare recessive neurometabolic disorder causing developmental delay, epilepsy, hypotonia, and dysmorphic features.
------	--

基因ID	728294
基因名	D2HGDH
Swiss	Q8N465
别名	D2HGDH;D2HGD

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn