

EIF6 Rabbit pAb

货号: B15134

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	23kDa/26kDa
实测分子量	27KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,293F,C2C12,Mouse liver,Rat thymus
细胞定位	Cytoplasm,Nucleus,nucleolus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-245 of human EIF6 (N P_852133.1).
序列	MAVRASFENNCEIGCFAKLTNTYCLVAIGGSENFYSVFEGELSDTIPVVHASIAGCRIIGRMCVGNRHHGLLVPNNTTDQELQ HIRNSLPDTVQIRRV EERLSALGNV TTCNDYVALVHPDL DRETEEILADV LKVEVFRQTVADQVLVGSYCVFSNQGGLVHP KTSIEDQDELSSLLQVPLVAGTVNRGSEVIAAGMVVNDWCAFCGLDTTSTELSVVESVFKLNEAQPSTIATSMRDSLIDSL T

靶点信息

研究背景	Hemidesmosomes are structures which link the basal lamina to the intermediate filament cytoskeleton. A n important functional component of hemidesmosomes is the integrin beta-4 subunit (ITGB4), a protein c ontaining two fibronectin type III domains. The protein encoded by this gene binds to the fibronectin type III domains of ITGB4 and may help link ITGB4 to the intermediate filament cytoskeleton. The encoded prot ein, which is insoluble and found both in the nucleus and in the cytoplasm, can function as a translation i nitiation factor and prevent the association of the 40S and 60S ribosomal subunits. Multiple non-protein c oding transcript variants and variants encoding two different isoforms have been found for this gene.
基因ID	3692
基因名	EIF6
Swiss	P56537
别名	EIF6;CAB;EIF3A;ITGB4BP;b(2)gcn;eIF-6;p27(BBP);p27BBP

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn