

RXRG Rabbit pAb

货号: B15136

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100
理论分子量	50kDa
实测分子量	61kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,DU145,U-87MG,Mouse brain,Rat brain
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-137 of human RXRG (NP_008848.1).
序列	MYGNYSHFMKFPAGYGGSPGHTGSTMSPSAALSTGKPMDSHPSTYDTPVSAPRTL SAVGTPLNALGSPYRVITSAMGP PSGALAAPPGINLVAPPSSQLNVVNSVSSSEDIKPLPLPGIGNMNYPSTSPGSLVKH

靶点信息

研究背景	This gene encodes a member of the retinoid X receptor (RXR) family of nuclear receptors which are involved in mediating the antiproliferative effects of retinoic acid (RA). This receptor forms dimers with the retinoic acid, thyroid hormone, and vitamin D receptors, increasing both DNA binding and transcriptional function on their respective response elements. This gene is expressed at significantly lower levels in non-small cell lung cancer cells. Alternatively spliced transcript variants have been described.
------	--

基因ID	6258
基因名	RXRG
Swiss	P48443
别名	RXRG;NR2B3;RXRC

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn