

Phospho-GAB1-Y627 Rabbit pAb

货号: B15146

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:100 - 1:200
理论分子量	76kDa/80kDa
实测分子量	125kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293
细胞定位	cell-cell junction,cytoplasm,cytosol
纯化	Affinity purification

抗原信息


抗原信息	A synthetic phosphorylated peptide around Y627 of human GAB1 (NP_002030.2?).
序列	VEYLD

靶点信息

研究背景	The protein encoded by this gene is a member of the IRS1-like multisubstrate docking protein family. It is an important mediator of branching tubulogenesis and plays a central role in cellular growth response, tr ansformation and apoptosis. Two transcript variants encoding different isoforms have been found for this gene.
基因ID	2549

基因名	GAB1
Swiss	Q13480
别名	GAB1;GRB2

产品验证



Western blot analysis of Phospho-GAB1-Y627 expressed in 293 using Phospho-GAB1-Y627 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn