

THAP2 Rabbit pAb

货号: B15163

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	26kDa
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse small intestine,Rat brain
细胞定位	nucleolus,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 83-228 of human THAP2 (NP_113623.1).
序列	CTHIKSMKLKSRNLLKKNNSCSPAGPSNLKSNISSQVLEHSYAFRNPMEAKKRIIKLEKEIASLRRKMKTCLQKERRATRR WIKATCLVKNLEANSVLPKGTSEHMLPTALSSLPLEDFKILEQDQDQDKTLLSLNLKQTKSTFI

靶点信息

研究背景	Predicted to enable DNA binding activity and metal ion binding activity. Located in nucleolus.
基因ID	83591
基因名	THAP2

Swiss	Q9H0W7
别名	THAP2

产品验证



Western blot analysis of THAP2 expressed in Mouse small intestine,Rat brain using THAP2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blotting buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn