

# RUVBL1/TIP49A/PONTIN Rabbit pAb

货号: B15208

产品信息

|       |                                                                                                             |
|-------|-------------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse                                                                                                 |
| 宿主    | Rabbit                                                                                                      |
| 克隆性   | Polyclonal                                                                                                  |
| 预测反应  |                                                                                                             |
| 应用    | WB IHC IF/ICC                                                                                               |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:200<br><b>IF/ICC:</b> 1:20 - 1:50                         |
| 理论分子量 | 42kDa/50kDa                                                                                                 |
| 实测分子量 | 55kDa                                                                                                       |
| 形式    | Liquid                                                                                                      |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.      |
| 偶联物   | Unconjugated                                                                                                |
| 阳性对照  | HL-60,MCF7,HepG2,293T,Mouse spleen,Mouse skeletal muscle,Mouse thymus,Mouse brain                           |
| 细胞定位  | Cytoplasm,Membrane,Nucleus,Nucleus matrix,centrosome,cytoskeleton,microtubule organizing center,nucleoplasm |
| 纯化    | Affinity purification                                                                                       |


抗原信息

|      |                                                                                                                                                                                                                                                |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human RUVBL1/TIP49A/PONTIN (NP_003698.1).                                                                                                               |
| 序列   | MKIEEVKSTTKTQRIASHSHVKGLGLDESLAKQAASGLVGQENAREACGVIVELIKSKKMAGRAVLLAGPPGTGKTALALAIAQELGSKVPFCPMVGSEVYSTEIKKTEVLMENFRRAIGLRIKETKEVYEGEVTETPCETENPMGGYGKTISHVIIGLKTAKGTKQLKLDPSIFESLQKERVEAGDVIYIEANS GAVKRQGRCDTYATEFDLEAEYVPLPKGDVHKKKEIIQDVTL |

靶点信息

|       |                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|-------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景  | This gene encodes a protein that has both DNA-dependent ATPase and DNA helicase activities and belongs to the ATPases associated with diverse cellular activities (AAA+) protein family. The encoded protein associates with several multisubunit transcriptional complexes and with protein complexes involved in both ATP-dependent remodeling and histone modification. Alternate splicing results in multiple transcript variants. |
| 基因ID  | 8607                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| 基因名   | RUVBL1                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Swiss | Q9Y265                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| 别名    | RUVBL1;ECP-54;ECP54;INO80H;NMP 238;NMP238;PONTIN;Pontin52;RVB1;TIH1;TIP49;TIP49A                                                                                                                                                                                                                                                                                                                                                       |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)