

# CELF2 Rabbit pAb

货号: B15217

## 产品信息

|       |                                                                                                        |
|-------|--------------------------------------------------------------------------------------------------------|
| 反应    | Human                                                                                                  |
| 宿主    | Rabbit                                                                                                 |
| 克隆性   | Polyclonal                                                                                             |
| 预测反应  |                                                                                                        |
| 应用    | WB                                                                                                     |
| 推荐浓度  | <b>WB:</b> 1:200 - 1:500                                                                               |
| 理论分子量 | 52kDa/54kDa/55kDa                                                                                      |
| 实测分子量 | 60kDa                                                                                                  |
| 形式    | Liquid                                                                                                 |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                           |
| 阳性对照  | MCF7,HeLa                                                                                              |
| 细胞定位  | Cytoplasm,Nucleus                                                                                      |
| 纯化    | Affinity purification                                                                                  |

## 抗原信息


|      |                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 422-521 of human CELF 2 (NP_006552.3). |
| 序列   | QQSAAGSQKEGPEGANLFYHLPQEFQDQDILQMFMFPGNVISAKVFIDKQTNLSKCFGVSYDNPVSAQAAIQAMNGFQIGMKRLKVQLKRSKNDSKPY                   |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | Members of the CELF/BRUNOL protein family contain two N-terminal RNA recognition motif (RRM) domain s, one C-terminal RRM domain, and a divergent segment of 160-230 aa between the second and third RR M domains. Members of this protein family regulate pre-mRNA alternative splicing and may also be involv ed in mRNA editing, and translation. Alternative splicing results in multiple transcript variants encoding di fferent isoforms. |
| 基因ID | 10659                                                                                                                                                                                                                                                                                                                                                                                                                                           |

|       |                                                      |
|-------|------------------------------------------------------|
| 基因名   | CELF2                                                |
| Swiss | O95319                                               |
| 别名    | CELF2;BRUNOL3;CELF-2;CUG-BP2;CUGBP2;ETR-3;ETR3;NAPOR |

产品验证



Western blot analysis of CELF2 expressed in MCF7,HeLa using CELF2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)