

ETV7 Rabbit pAb

货号: B15223

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	22kDa/29kDa/30kDa/32kDa/36kDa/38kDa
实测分子量	37kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-937,A-431,U-251MG,HeLa
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human ETV7 (NP_057219.1).
序列	MQEGELAISPISPAAMPPLGTHVQARCEAQINLLGEGGICKLPGRRLRIQPALWSREDVLHWLRWAEQEYSLPCTAEHGFE MNGRALCILT KDDFRHRAPSSGDVLYELLQYIKTQRRALVCGPFFGGIFRLKTPTQHSPVPPEEVTGPSQMDTRRGHLLQP PDPGLTSNFGHLDDPGLARWTPGKEESLNLCHCAELGCRTQGVCSPFA

靶点信息

研究背景	The protein encoded by this gene belongs to the ETS family of transcription factors, which is a large group of evolutionarily conserved transcriptional regulators that play an important role in a variety of cellular processes throughout development and differentiation, and are involved in oncogenesis as well. This protein is predominantly expressed in hematopoietic tissues. Several alternatively spliced transcript variants encoding different isoforms have been described for this gene (PMID:11108721).
------	---

基因ID	51513
基因名	ETV7
Swiss	Q9Y603
别名	ETV7;TEL-2;TEL2;TELB

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn