

# PCGF1 Rabbit pAb

货号: B15227

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                      |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB IF/ICC                                                                                            |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:1000<br><b>IF/ICC:</b> 1:50 - 1:200                                             |
| 理论分子量 | 20kDa/30kDa                                                                                          |
| 实测分子量 | 35KDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | HeLa,A-549,MCF7,Mouse testis,Mouse brain,Rat testis                                                  |
| 细胞定位  | Nucleus                                                                                              |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 160-259 of human PCGF 1 (NP_116062.2). |
| 序列   | AHYRYDEQLNLCLERLSSGKDKNKSVLQNKYVRCSVRAEVRHLRRVLCHRLMLNPQH VQLLF DNEVL PDHMTMKQIWL SRWFGKPSPLLLQYSVKEKRR              |

靶点信息

|      |                                                                                                                                                                                                                                                                                        |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | PCGF1 is a mammalian homolog of the Drosophila polycomb group genes, which act as transcriptional repressors to regulate anterior-posterior patterning in early embryonic development (Nunes et al., 2001 [PubMed 11287196]). See also PCGF2 (MIM 600346).[supplied by OMIM, Aug 2008] |
| 基因ID | 84759                                                                                                                                                                                                                                                                                  |

|       |                                         |
|-------|-----------------------------------------|
| 基因名   | PCGF1                                   |
| Swiss | Q9BSM1                                  |
| 别名    | PCGF1;2010002K04Rik;NSPC1;RNF3A-2;RNF68 |

产品验证



Western blot analysis of PCGF1 expressed in Mouse testis,Rat testis,HeLa using PCGF1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)