

KYNU Rabbit pAb

货号: B15243

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:100 - 1:500 IHC: 1:50 - 1:100 IF/ICC: 1:50 - 1:100
理论分子量	34kDa/52kDa
实测分子量	52KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,HAP1(Negative control)
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-307 of human KYNU (NP_001028170.1).
序列	MEPSSLELPADTVQRIAAELKCHPTDERVALHLDEEDKLRHFRFCFYIPKIQDLPPVDLSLVNKDENAIYFLGNSLGLQPKMV KTYLEEEELDKWAKIAAYGHEVGKRPWITGDESIVGLMKDIVGANEKEIALMNALTVDNLHLLMSFFKPTPKRYKILLEAKAFP SDHYAIESQLQLHGLNIEESMRMIKPREGEETLRIEDILEVIEKEGDSIAVILFSGVHFYTGQHFNIPAITKAGQAKGCYVGFDL AHAVGNVELYLHDWGVDFACWCSYKYLNAGAGGIAGAFIHEKHAHTIKPARSEFFN

靶点信息

研究背景	Kynureninase is a pyridoxal-5'-phosphate (pyridoxal-P) dependent enzyme that catalyzes the cleavage of L-kynurenine and L-3-hydroxykynurenine into anthranilic and 3-hydroxyanthranilic acids, respectively. Ky nureninase is involved in the biosynthesis of NAD cofactors from tryptophan through the kynurenine path way. Alternative splicing results in multiple transcript variants.
------	--

基因ID	8942
基因名	KYNU
Swiss	Q16719
别名	KYNU;KYNUU

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn