

ASMTL Rabbit pAb

货号: **B15269**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IP
推荐浓度	WB: 1:500 - 1:2000 IP: 1:50 - 1:200
理论分子量	62kDa/67kDa/68kDa
实测分子量	72kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,22Rv1,SW620,NIH/3T3,HeLa,Mouse pancreas,Mouse thymus,Mouse testis
细胞定位	cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human ASMTL (NP_004183.2).
序列	MVLCPVIGKLLHKRVVLASASPRRQEILSNAGLRFEVVPSEKFKELDKASFATPYGYAMETAKQKALEVANRLYQKDLRAPD VVIGADTIVTVGGLILEKPVDKQDAYRMLSRLSGREHSVFTGVAIVHCSSKDHQLDTRVSEFYEETKVKFSEELSELLWEYV HSGEPMDKAGGYGIQALGGMLVESVHGDFLNVVGFLNHFCKQLVKLYPPRPEDLRRSVKHDSIPAADTFEDLSDVEG GGSEPTQRDAGSRDEKAEAGEAGQATAEAECHRTRETLPFPFTRLLE

靶点信息

研究背景	The protein encoded by this gene has an N-terminus that is similar to the multicopy associated filamentat ion (maf) protein of Bacillus subtilis and to orfE of Escherichia coli, while the C-terminus is similar to N-acetylserotonin O-methyltransferase. This gene is located in the pseudoautosomal region 1 (PAR1) of X and Y chromosomes. Three transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	8623
基因名	ASMTL
Swiss	O95671
别名	ASMTL;ASMTLX;ASMTLY;ASTML

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn