

LHX6 Rabbit pAb

货号: B15285

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	19kDa/38kDa/40kDa/41kDa/43kDa
实测分子量	45kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW480,U-87MG,HepG2,NIH/3T3,Mouse brain,Mouse liver,Mouse heart,Rat brain,Rat heart
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 143-392 of human LHX6 (NP_055183.2).
序列	NKEIFCKMDYFSRFGTKCARCGRQIYASDWVRRARGNAYHLACFACFSCKRQLSTGEEFGLVEEKVLCRIHYDTMIENLKR AAENGNGLTLEGAVPSEQDSQPKPAKRARTSFTAEQLQVMAQFAQDNNPDAQTLQKLADMTGLSRRVIQVWFQNCR ARHKKHTPQHVPVPSGAPPSRLPSALSDDIHYTPFSSPERARMVTLHGYIESQVQCGQVHCRLPYTAPPVHLKADMDGPL SNRGEKVILFQY

靶点信息

研究背景	This gene encodes a member of a large protein family that contains the LIM domain, a unique cysteine-rich zinc-binding domain. The encoded protein has tandem LIM domains as well as a DNA-binding homeodomain. The protein functions as a transcription factor involved in embryogenesis and head development and is highly expressed in neural crest derived mesenchyme cells. Alternative splicing results in multiple transcript variants encoding distinct isoforms.
------	---

基因ID	26468
基因名	LHX6
Swiss	Q9UPM6
别名	LHX6;LHX6.1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn