

HLX Rabbit pAb

货号: B15289

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	50kDa
实测分子量	51kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HL-60,Raji,THP-1,293T
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 80-280 of human HLX (NP_068777.1).
序列	HPHASFQAAARSPLRPTPVVAPSEVPAGFPQRLSPLSAAYHHHPQQQQQQPQQQPPPPPRAGALQPPASGTRVV PNPHHSGSAPAPSSKDLKFGIDRILSAEFDPKVKEGNTLRDLTSLTGGRPAGVHLSGLQPSAGQFFASLDPINEASAILSP LNSNPRNSVQHQQFQDTFPGPYAVLTKDTMPQTYKRKRWSR

靶点信息

研究背景	Transcription factor required for TBX21/T-bet-dependent maturation of Th1 cells as well as maintenance of Th1-specific gene expression. Involved in embryogenesis and hematopoiesis (By similarity.
基因ID	3142

基因名	HLX
Swiss	Q14774
别名	HLX;HB24;HLX1

产品验证



Western blot analysis of HLX expressed in HL-60,Raji,THP-1 using HLX Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn