

VAMP3/Cellubrevin Rabbit pAb

货号: B15290

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	11kDa
实测分子量	15kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SKOV3,U-251MG,SW480,A-549,293T,Mouse testis,Mouse kidney,Mouse lung,Rat brain
细胞定位	Cell junction,Membrane,Single-pass type IV membrane protein,synapse,synaptosome
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human VAMP3/C ellubrevin (NP_004772.1).
序列	MSTGPTAATGSNRRLLQQTQNQVDEVVDIMRVNVDKVLERDQKLSEDDRADALQAGASQFETSAAKLKRKYWWKNCK MWA

靶点信息

研究背景	Synaptobrevins/VAMPs, syntaxins, and the 25-kD synaptosomal-associated protein are the main components of a protein complex involved in the docking and/or fusion of synaptic vesicles with the presynaptic membrane. This gene is a member of the vesicle-associated membrane protein (VAMP)/synaptobrevin family. Because of its high homology to other known VAMPs, its broad tissue distribution, and its subcellular localization, the protein encoded by this gene was shown to be the human equivalent of the rodent cellubrevin. In platelets the protein resides on a compartment that is not mobilized to the plasma membrane on calcium or thrombin stimulation.
------	---

基因ID	9341
基因名	VAMP3
Swiss	Q15836
别名	VAMP3;CEB

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn