

ZNF346 Rabbit pAb

货号: B15293

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	29kDa/32kDa/35kDa
实测分子量	33-38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,MCF7,THP-1,HeLa,Mouse testis,Mouse heart,Mouse brain
细胞定位	Cytoplasm,Nucleus,nucleolus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-294 of human ZNF346 (NP_036411.1).
序列	MEYPAPATVQAADGGAAGPYSSSELLEGQEPDGVRFDREARRRLWEAVSGAQPVGREEVEHMIQKNQCLFTNTQCKVC CALLISESQKLAHYQSKKHANKVKRYLAIHGMETLKGETKKLSDQKSSRSKDKNQCCPICNMTFSSPVVAQSHYLGKTH AKNLKQKQSTKVEALHQNREIDPDKFCSLCHATFNDPVMAQQHYVGKKHRKQETKLKLMARYGRLADPAVTDFPAGK GYPCKTCKIVLNSIEQYQAHVSGFKHNQSPKTVASSLGQIPMQRQPIQKDSTTLED

靶点信息

研究背景	The protein encoded by this gene is a nucleolar, zinc finger protein that preferentially binds to double-stranded (ds) RNA or RNA/DNA hybrids, rather than DNA alone. Mutational studies indicate that the zinc finger domains are not only essential for dsRNA binding, but are also required for its nucleolar localization. The encoded protein may be involved in cell growth and survival. It plays a role in protecting neurons by inhibiting cell cycle re-entry via stimulation of p21 gene expression. Alternative splicing of this gene results in multiple transcript variants.
------	--

基因ID	23567
基因名	ZNF346
Swiss	Q9UL40
别名	ZNF346;JAZ;Zfp346

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn