

BOLL Rabbit pAb

货号: B15316

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	19kDa/31kDa/32kDa/37kDa
实测分子量	35-42KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse testis,Rat testis
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-283 of human BOLL (NP_149019.1).
序列	MQTDSLSPSPNPVSPVPLNNPTSAPRYGTVIPNRIFVGGIDFKTNESDLRKFFSQYGSVKEVKIVNDRAGVSKGYGFVTFET QEDAQKILQEAELNYKDKKLNIGPAIRKQQVGI PRSSIMPAAGTMYLTTSTGYPTYHNGVAYFHTPEVTSVPPPWP SRSV CSSPVMVAQPIYQQPAYHYQATTQYLPGQWQWSVPQPSASSAPFLYLQPSEVIYQPVEIAQDGGCVPPPLSLMETS VPEP YSDHGVQATYHQVYAPSAITMPAPVMQPEPIKTVWSIH Y

靶点信息

研究背景	This gene belongs to the DAZ gene family required for germ cell development. It encodes an RNA-binding protein which is more similar to Drosophila Boule than to human proteins encoded by genes DAZ (deleted in azoospermia) or DAZL (deleted in azoospermia-like). Loss of this gene function results in the absence o f sperm in semen (azoospermia). Histological studies demonstrated that the primary defect is at the meio tic G2/M transition. Two alternatively spliced transcript variants encoding distinct isoforms have been fou nd for this gene.
------	---

基因ID	66037
基因名	BOLL
Swiss	Q8N9W6
别名	BOLL;BOULE

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn