

PPP4R1 Rabbit pAb

货号: B15319

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IP
推荐浓度	WB: 1:500 - 1:2000 IP: 1:50 - 1:100
理论分子量	105kDa/107kDa
实测分子量	130kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,U-251MG,BT-474,K-562,SKOV3
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 324-647 of human PPP4 R1 (NP_001035847.1).
序列	NPSSSGQYFKEESKSSEEMSVENKNRTRDQEAPEDVQVRPEDTPSDLSVSNSSVILENTMEDHAAEASGKPLGEISVPLD SSLLCTLSSESHQEAASNENDKKPGNYKSMLRPEVGTTSQDSALLDQELYNSFHFWRTPLPEIDLDIELEQNSGGKPSPEG PEEESEGPVPSSPNITMATRKELEEMIENTLEPHIDDPDVKAQVEVLSAALRASSLDAHEETISIEKRSDLQDELDINELPNCKI NQEDSVPLISDAVENMDSTLHYIHSDSDLSNNSFSFSPDEERRTKVQDVVPQALLDQYLSMTDPSRAQTVDTETIAKHCA

靶点信息

研究背景	This gene encodes one of several alternate regulatory subunits of serine/threonine protein phosphatase 4 (PP4). The protein features multiple HEAT repeats. This protein forms a complex with PP4RC. This complex may have a distinct role from other PP4 complexes, including regulation of HDAC3 (Zhang et al., PMID: 15805470). There is also a transcribed pseudogene on chromosome 20. Alternative splicing results in multiple transcript variants.
------	---

基因ID	9989
基因名	PPP4R1
Swiss	Q8TF05
别名	PPP4R1;MEG1;PP4(Rmeg);PP4R1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn