

CRHR1 Rabbit pAb

货号: B15321

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	28kDa/43kDa/46kDa/47kDa/50kDa
实测分子量	60kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,SH-SY5Y,HT-29,HeLa,Mouse brain,Mouse liver,Mouse testis,Rat brain
细胞定位	Cell membrane,Endosome,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 24-120 of human CRHR 1 (NP_001138618.1).
序列	SLQDQHCESLSLASNISGLQCNASVDLIGTCWPRSPAGQLVVRPCPAFFYGVRynttNNGYRECLANGSWAARVNYSECQEILNEEKKSkvHYHVAV

靶点信息

研究背景	This gene encodes a G-protein coupled receptor that binds neuropeptides of the corticotropin releasing hormone family that are major regulators of the hypothalamic-pituitary-adrenal pathway. The encoded protein is essential for the activation of signal transduction pathways that regulate diverse physiological processes including stress, reproduction, immune response and obesity. Alternative splicing results in multiple transcript variants. Naturally-occurring readthrough transcription between this gene and upstream GeneID:147081 results in transcripts that encode isoforms that share similarity with the products of this gene.
------	--

基因ID	1394
基因名	
Swiss	P34998
别名	CRHR1;CRF-R;CRF-R-1;CRF-R1;CRF1;CRFR-1;CRFR1;CRH-R-1;CRH-R1;CRHR;CRHR1L

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn