

MYO10 Rabbit pAb

货号: B15323

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	11kDa/163kDa/237kDa
实测分子量	270kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,Mouse lung
细胞定位	Cell projection,Cytoplasm,Peripheral membrane protein,cell cortex,cytoskeleton,cytosol,filopodium membrane,filopodium tip,lamellipodium,ruffle
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 845-944 of human MYO10 (NP_036466.2).
序列	EAELRAQQEEETRKQQELEALQKSQKEALTRELEKQKENKQVEEILRLEKEIEDLQRMKEQQELSLTEASLQKLQERRDQELRRLEEEACRAAQEFLES

靶点信息

研究背景	This gene encodes a member of the myosin superfamily. The protein represents an unconventional myosin; it should not be confused with the conventional non-muscle myosin-10 (MYH10). Unconventional myosins contain the basic domains of conventional myosins and are further distinguished from class members by their tail domains. This gene functions as an actin-based molecular motor and plays a role in integration of F-actin and microtubule cytoskeletons during meiosis.
------	--

基因ID	4651
基因名	MYO10
Swiss	Q9HD67
别名	MYO10;myosin X

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn