

ERVFRD-1 Rabbit pAb

货号: **B15340**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	59kDa
实测分子量	60KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Cell membrane,Single-pass membrane protein,Virion
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 16-250 of human ERVFR D-1 (NP_997465.1).
序列	AYRHPDFPLLEKAQQLQSTGSPYSTNCWLCTSSSTETPGTAYPASPREWTSIEAELHISYRWDPNLKGLMRPANSLLSTV KQDFPDIRQKPPIFGPIFTNINLMGIAPICVMAKRKNGTNVGTLPSTVCNVTFFTVDSNQQTYYTYTHNQFRHQPRFPKPPNI TFPQGTLLDKSSRFCQGRPSSCSTRNFWFRPADYNQCLQISNLSSTAIEWVLLDQTRNSLFWENKTKGANQSQ

靶点信息

研究背景	Many different human endogenous retrovirus (HERV) families are expressed in normal placental tissue at high levels, suggesting that HERVs are functionally important in reproduction. This gene is part of a human endogenous retrovirus provirus on chromosome 6 that has inactivating mutations in the gag and pol genes. This gene is the envelope glycoprotein gene which appears to have been selectively preserved. The gene's protein product plays a major role in placental development and trophoblast fusion. The protein has the characteristics of a typical retroviral envelope protein, including a cleavage site that separates the surface (SU) and transmembrane (TM) proteins which form a heterodimer.
基因ID	405754
基因名	ERVFRD-1
Swiss	P60508
别名	ERVFRD-1;ERVFRDE1;GLLL6191;HERV-FRD;HERV-W/FRD;UNQ6191;envFRD

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn